

phantom Limb

Tech Rider




Phantom Limb (2023) is a project exploring the enigmatic and poetic relationship between a human being and the black box that is their interior through the use of a medical ultrasound machine. There are two versions of the work:

**The performance**

Uses the ultrasound as a source for live video and sound and least around 15 minutes, the performance includes the artist, ultrasound machine, PA, projector, lights and Blue curtains.

**The installation**


Is a stand-alone, using the ultrasound machine to paly back video through CRT and sound through headphones.




Contact

Amos Peled  
amos.peled97@gmail.com  
+31623371817  
amospeled.net

# Performance set up



- Red = Provided by host
- Green = Provided by artist




## Performance

Name	Details	Amount	RED - HOST GREEN - ARTIST
<b>Vidoe and lighth</b>			
BLACKMAGIC ATEM MINI PRO	HDMI to USB	1	iii
Blackmagic intensity shutter	HDMI to USB	1	
Extron IN1604 HD	VGA to HDMI	1	
Back projector short throw	4000 lumens	1	
Projection screen	4:3 aspect ratio	1	
Theater light with mounting	<b>To be discussed</b>	1	
Spaces Lighths	LED DMX control	1	
DMX control unit	Small mixer for room lighth	1	
<b>Vidoe and lighth cables</b>			
Details		Amount	
HDMI	10 M	1	
HDMI	0.5 M	2	
VGA	0.5 M	1	
DMX	10 M	4	
DMX	1 M	1	
<b>Audio</b>			
Details		Amount	
Audio interface		1	
PA speaker active		2	
Sub		1	
Stands for speakers		2	
DI	Stereo passive	1	
<b>Audio cables</b>			
Details		Amount	
Jack stereo to 2 mono jack			
XLR	10 M	1	
XLR	20 M	1	
<b>General</b>			
Details		Amount	
Ultrasound Machines		1	
Chairs	<b>To be discussed</b>	20-40	
Electricity	<b>As shown in the picture</b>		

**Note about electricity: min 2 different phases light and sound+set up preferences is 3 phase 1 light 2 sound 3 set up. The ultrasound is running on 2KV**

**The walls will be coverd by blue curtains that will be provided by the host, the artist is waiting for measurements of the space to decide on the curtains**

Picture of the performance



AP  
0.038

01 Feb 23

20:30:13

8L5 3Hz

8.0MHz 45mm

PHANTOM LIMB

Low Flow /V

S2/-3/ 5/E:1

3/3 00.40MHz

CD Gain = 1

CDE 26dB

RICA-DOP


Angle L/R

CDV/VE


CD Pwr

CD Size

# Installation set up



- Red = Provided by host
- Green = Provided by artist



## Installation

Name	Details	Amount	Link	Status
<b>Vidoe and Audio</b>				
CRT screen		1		Green
Headphones splitter amplifier		1		Red
Headphones		2		Red
Vidoe Player	S-video out, loop playback	1		Red
<b>General</b>				
Ultrasound Machines		1		Green
Electricity		1		Red

**Note about the installation:** the installation is intended to run for around a month, so the better to use a good quality video player that supports auto paly back after turning off, in case problem with the electricity



Picture of the installation